

Pipet Tips

I. APPLICATIONS

1. For use with pipettors for the measure and transfer of liquids/chemicals

II. FEATURES

1. Materials: Strong, medical-grade polypropylene
2. Chemically resistant to alcohols and mild organic solvents
3. 100% inert
4. Both tips and racks are autoclavable
5. Universal fit
6. Certified Rnase & Dnase Free

III. TIP DESIGN

1. 10µl, 200µl, and 1000µl tip sizes
2. All sizes available in non-sterile bulk or sterile racked packaging
3. 10µl and 200µl tip are graduated, 1000µl is non-graduated

IV. PACKAGING

1. Bulk Tips (10µl, 200µl, 1000µl, and Extended 1000µl)
 - a). Packaged 1,000 Tips/bag and 10 bags/case
2. Racked Tips (10µl, 200µl, 1000µl, and Extended 1000µl)
 - a). 10µl, 200µl: Packaged 96 tips/rack and 100 racks/case
 - b). 1000µl and Extended 1000µl: Packaged 100 tips/rack and 60 racks/case

V. DESCRIPTION

Item	Size	Packaging	Sterile
20-0010	10µl	Bulk, 1,000 tips/bag, 10 bags/case	No
20-0200	200µl	Bulk, 1,000 tips/bag, 10 bags/case	No
20-1000	1000µl	Bulk, 1,000 tips/bag, 10 bags/case	No
20-1250	Extended 1000µl	Bulk, 1,000 tips/bag, 10 bags/case	No
21-0010	10µl	Racked, 96 tips/rack and 100 racks/case	Yes
21-0200	200µl	Racked, 96 tips/rack and 100 racks/case	Yes
21-1000	1000µl	Racked, 100 tips/rack and 60 racks/case	Yes
21-1250	Extended 1000µl	Racked, 100 tips/rack and 60 racks/case	Yes

VI. STERILITY

1. Products are EO Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for

product release.

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


IX. PRODUCT PARAMETER

Dimensions = mm


Cat #	Description	A	B	C	D	E	F	G	H
20-0010	Pipet Tip	31.6	O.D. 6.0	O.D. 3.8	O.D. 0.86	3.9	16.8	24.7	-
20-0200	Pipet Tip	50.5	O.D. 7.5	O.D. 6.1	O.D. 0.9	16.1	-	-	-
20-1000	Pipet Tip	77.9	O.D. 9.6	O.D. 8.2	O.D. 1.1	20.3	-	-	-
20-1250	Pipet Tip	88.42	43.41	17.75	25.53	O.D. 7.43	O.D. 8.6	O.D. 1.3	O.D.0.72




20-0010




20-0200



20-1000





20-1250



21-0010

Dimensions=mm



Cat #	Description	A	B	C	D	E	F	G
21-0010	Rack	85.79	81.79	80.8	16.7	21.9	83.79	85.79
		H	I	J	K	L	M	-
		15.7	119.6	8.98	O.D.4.85	85.79	120.6	-



21-0200

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
21-0200	Rack	82	81.2	77.8	44.3	32.6	31.6	84
		H	I	J	K	L	M	-
		84	82	O.D.4.8	8.8	118	117.8	-



21-1000

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
21-1000	Rack	109.3	103.1	101.4	48.43	65.2	111.3	111.3
		H	I	J	K	L	M	-
		11.52	19.95	∅ 8.02	136	128	125.8	-

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
21-1250	Rack	69.8	30.2	26	3.2	18.6	191.91	173.8
		H	I	J	K	L	-	-
		85.2	85.4	99	63	9	-	-

