

Filter Tips

I. APPLICATIONS

1. For use with pipettors for the measure and transfer of liquids/chemicals

II. FEATURES

1. Materials: Tips are made of PP, filters are made of UHMWPE
2. Chemically resistant to alcohols and mild organic solvents
3. 100% inert
4. Both tips and racks are both autoclavable
5. Certified Rnase & Dnase Free

III. FILTER TIPS DESIGN

1. 10 μ l, 200 μ l, and 1000 μ l tip sizes
2. All sizes available in non-sterile bulk or sterile racked packaging
3. The filter's average pore size of 20 μ m eliminates errors caused by saturation vapor pressure

IV. PACKAGING

1. Bulk Tips
 - a). Packaged 1,000 Pieces/Pack and 10 Packs/case
2. Racked Tips
 - a). Packaged 96 Pieces/Rack, 10 Racks/Pack, 10 Packs/Case(10 μ l, 200 μ l,)
 - b). Packaged 96 Pieces/Rack, 6 Racks/Pack, 10 Packs/Case(1000 μ l)

V. DESCRIPTION

Item	Size	Packaging	Sterile
22-0010	10 μ L	1000 Pieces/Pack, 10 Packs/Case	No
22-0011	10 μ L	1000 Pieces/Pack, 10 Packs/Case	No
23-0010	10 μ L	96 Pieces/Rack, 10 Racks/Pack, 10 Packs/Case	Yes
23-0011	10 μ L	96 Pieces/Rack, 10 Racks/Pack, 10 Packs/Case	Yes
22-0200	200 μ L	1000 Pieces/Pack, 10 Packs/Case	No
22-0201	200 μ L	1000 Pieces/Pack, 10 Packs/Case	No
23-0200	200 μ L	96 Pieces/Rack, 10 Racks/Pack, 10 Packs/Case	Yes
23-0201	200 μ L	96 Pieces/Rack, 10 Racks/Pack, 10 Packs/Case	Yes
22-1000	1000 μ L	1000 Pieces/Pack, 10 Packs/Case	No
22-1001	1000 μ L	1000 Pieces/Pack, 10 Packs/Case	No
23-1000	1000 μ L	96 Pieces/Rack, 6 Racks/Pack, 10 Packs/Case	Yes
23-1001	1000 μ L	96 Pieces/Rack, 6 Racks/Pack, 10 Packs/Case	Yes

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER


Dimensions=mm

Cat #	Description	A	B	C	D	E	F
22-0010	Filter Tips	39.10	18.00	6.30	5.00	0.40	0.70
22-0011	Filter Tips	39.10	18.00	6.30	5.00	0.40	0.70




10µL Filter Tips

Cat #	Description	A	B	C	D	E	F
22-0200	Filter Tips	59.30	16.00	7.10	5.23	0.52	0.92
22-0201	Filter Tips	59.30	16.00	7.10	5.23	0.52	0.92





200µL Filter Tips



Cat #	Description	A	B	C	D	E	F
22-1000	Filter Tips	88.10	9.00	8.50	7.20	0.82	1.2
22-1001	Filter Tips	88.10	9.00	8.50	7.20	0.82	1.2



1000µL Filter Tips



23-0010 & 23-0011



23-0200 & 23-0201

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
23-1000	Rack	69.8	30.2	26	3.2	18.6	191.91	173.8
23-1001		H	I	J	K	L	-	-
		85.2	85.4	99	63	9	-	-

