

## Polypropylene Cryogenic Storage Boxes

### I. APPLICATIONS

1. Storage and organization of any range of 0.5ml to 2.0ml microtubes/cryogenic vials

### II. FEATURES

1. Materials: durable polypropylene
2. Chemically resistant to alcohols and mild organic solvents
3. Temperature range: stable from -80°C to 121°C (Autoclavable)

### III. BOX DESIGN

1. 9x9 well storage (81-wells total) or 10x10 well storage (100-wells total) or 5x5 well storage (25-wells total)
2. 5 assorted colors: clear, pink, green, blue and yellow
3. Holds 0.5ml, 1.5ml or 2.0ml microtubes
4. 81-well box features a removable lid includes a numbered grid for orientation
5. 90-9100 box features a 3-hinged lid with alphanumeric grid for orientation
6. 25-well box features a removable lid includes a numbered grid for orientation

### IV. PACKAGING

1. 81-well and 100-well: Bulk, 20 boxes/case, assorted colors  
25-well: Bulk, 30 boxes/case, assorted colors

### V. DESCRIPTION

Item	Description	Packaging	Sterile
90-9081	Polypropylene 81-well Cryogenic storage box	Bulk, 5 /bag, 4 bags/case, assorted colors	No
90-9100	Polypropylene 100-well Cryogenic storage box with a 3-hinged lid	Bulk, 5 /bag, 4 bags/case, assorted colors	No
90-9200	Polypropylene 100-well Cryogenic storage box	Bulk, 5 /bag, 4 bags/case, assorted colors	No
90-9025	Polypropylene 25-well Cryogenic storage box	Bulk, 5 /bag, 6 bags/case, assorted colors	No

### VI. STERILITY

1. N/A

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


## VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.


## IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-9081	Box	132.0	12.2	35.0	30.0	26.0	16.5	132.0
90-9081C	Lid	132.0	132.0	26.0	-	-	-	-






90-9081



90-9081C




Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G	H	I
90-9100	Box	140	12.2	140	12.2	54.0	50.0	32.0	90.0	24.0






**Dimensions=mm**

Cat #	Description	A	B	C	D	E	F	G
90-9200	Box	132	129.7	130.7	12.6	31.3	26.4	4.9
		H	I	J	K	L		
90-9200C	Lid	132.5	26.7	25.0	38.4	84.7		




Dimensions=mm

Cat #	Description	A	B	C	D	E	F	G
90-9025	Box	75.7	12.2	35.0	30.0	26.0	16.5	75.7
90-9025C	Lid	75.7	75.7	26.0	-	-	-	-



90-9025



90-9025C