

## Silicone Sealing Mat 02-1096

### I. APPLICATIONS

1. Made of silicon, having a temperature resistance of at least -80 ° C to + 50 ° C
2. To seal square deep well plates, 96 wells, with 96 small protrusions integrated
3. Cannot be punctured

### II. PACKAGING

1. Bulk, non-sterile.
2. 10 PCS/Bag, 5 Bags/Case

### III. CERTIFICATION


Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

### VI. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

### V. PRODUCT PARAMETER

Dimensions = mm



Item	Product dimension				Well dimension		Theoretical Weight
	Length	Width	Height	Thickness	Edge length	Position	
	L	W	H	t	A/B	P1/P2	
02-1096	123.00	81.00	6.00	1.00	8.30	9.00	31.9